

NOW

PRESSURE
RATED AND
ULTRASOUND
READY

More tools for the toolbox

The
M/29[®] midterm[®]

Midline Catheter System



The Right Line from the Start


FLEXICATH

More tools for the toolbox

The M/29[®] midterm[®] Midline Catheter System

NOW
PRESSURE
RATED AND
ULTRASOUND
READY



The M/29 Midterm[®] is the world's first pressure injectable IV catheter solution to provide a "built in" 360° Maximum Barrier[®] against airborne and touch borne contamination during the insertion process and offers a simple to place intermediate dwell alternative for patients needing IV therapy for up to 30 days. The M/29 product suite now includes ultrasound compatible, pressure injectable devices (see table below) that extend the application of the M/29 in the clinical setting to include ultrasound assisted placement and contrast agent administration – providing additional tools for the vascular access specialist, improving patient satisfaction, and enabling better care.

The M29[®] was designed to provide a safe, easy-to-insert pressure injectable catheter for intermediate dwell applications. The M29[®] is the only touch-free midline that can be placed in under 10 minutes following minimal training.

Product Code	Description	Catheter Gauge	Catheter Length	Needle Length	Flow Rate	Packaging
3SP0101	Kit - M/29 SNM Pressure Injectable Catheter System	20	19 CM	5 CM	5 mL/sec	5 per Carton
3SP1101	Kit - M/29 SNM Pressure Injectable Catheter System	20	19 CM	7 CM	5 mL/sec	5 per Carton
3SP0103	Kit - M/29 SNM Pressure Injectable Catheter System	20	10 CM	5 CM	5 mL/sec	5 per Carton
3SP1103	Kit - M/29 SNM Pressure Injectable Catheter System	20	10 CM	7 CM	5 mL/sec	5 per Carton
3N0102	Introducer, OTN, 19 Gauge, 5CM	-	-	5 CM	-	5 per Carton
3N1102	Introducer, OTN, 19 Gauge, 5CM	-	-	7 CM	-	5 per Carton

Therapies not appropriate for the M/29 Midterm[®] catheter include: continuous vesicants, parenteral nutrition, solutions with pH<5 and >9, solutions with osmolarity greater than 600 mOsm/L.

FLEXICATH

www.flexicath.com

Dealer imprint:

M/29, Midterm, PeelGuard, 3D Maximum Barrier, and 360 Maximum Barrier are registered trademarks of Flexicath.
© 2011, Flexicath. All rights reserved.



Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

MKT139 rev.1.0
October 08, 2012